



# **MicroScale Conductivity Temperature Instrument**

Model 125

**Operator's Manual**  
(Revision B 02-DEC-2008)

# MicroScale Conductivity-Temperature Instrument Model 125

## CONTENTS

INTRODUCTION  
Dos and DON'Ts  
SPECIFICATIONS  
PANEL FUNCTIONS

GROUNDING  
CALIBRATION  
MAINTENANCE  
REFERENCES

## INTRODUCTION

The Model 125 MicroScale Conductivity-Temperature Instrument (MSCTI) is designed to measure the temperature and electrical conductivity of water solutions containing conductive ions. The MSCTI provides analog voltage outputs that are functions of the solution's electrical conductivity and temperature. The instrument is intended for use in moving solutions where spatial resolution and time response are important.

## DOs AND DON'Ts

Always use fresh, clean solutions. Large contaminants such as algae flakes can become caught on the sensor's electrodes and cause calibration shifts. If possible, position the sensor so that it points directly upstream in the flow. This orientation is less susceptible to contamination.

Debris caught on the electrodes can sometimes be removed by gently rinsing them under a faucet or using a squeeze bottle. Never tap the stainless steel shaft. This causes vibrations that can cause the glass sensor to become cracked.

Avoid solutions that contain oil or slime. The conductivity electrodes are platinized to increase their ability to pass electrical current into the solution. This plating is sponge-like. If the plating becomes fouled with oil or slime, then the electrodes will have difficulty conducting electricity. The 'PHASE ERROR' light on the electronic bridge's front panel comes on when this condition occurs.

Minimize the time the sensor is exposed to air after calibration. Calibration shifts will occur if the air exposure time is longer than a few seconds. When moving the sensor from one solution to another, take advantage of the fact that a small drop of solution will cling to the sensor's electrodes and protect them.

Always rinse the sensor with fresh water after use.

Never expose the sensor to water unnecessarily.

The preamplifier is watertight and may be submerged. The cable may also be submerged. The connector to the electronic bridge is not watertight. Do not allow water to get into this connector. If this happens, then the instrument will not operate properly and can be permanently damaged. When storing or transporting a wet sensor assembly, be sure that the connector is kept at the highest point. Small drops of water tend to run down along the cable and into the connector if allowed.

The preamplifier cable is a complicated instrumentation cable. It is important that the cable jacket not allow water to enter the cable. Never allow the cable to be damaged. Don't allow it to lay on the floor where it can be stepped on or have equipment rolled over it.

## SPECIFICATIONS

CONDUCTIVITY	
Measurement equation	$V_o = G * C + V_{off}$
8 hour calibration stability	Better than 1% of C reading
Time response	-3db at approx. 800 hz
Spatial response	-3db at approx. 4 cyc/cm
Noise	< 1mV RMS 10 Hz to 1 kHz
Output voltage range	+/- 5 volts
Linear conductivity range	0.5 to 800 mS/cm

TEMPERATURE	
Measurement equation	$V_o = G \exp (A+B/T) + V_{off}$
8 hour calibration stability	Better than 0.01 deg C
Time response	Approx. 7E-3 seconds
Noise	< 1mV RMS 10 Hz to 200 Hz
Output voltage range	+/- 5 volts
Temperature range	-10 to 100 deg C

## PANEL FUNCTIONS

The Model 125 front panel has several displays and controls. Each of these is described below:

**POWER** - This switch controls power to the electronic bridge and preamplifier assembly. The adjacent light is 'on' when power is applied. Note that the switch simply interrupts power from the power supply module and does not control the power supply module in any way.

**GAIN - TEMPERATURE** - This controls 'G' in the temperature measurement equation mentioned above. 'CW' operation of this adjustment increases 'G'.

**GAIN - CONDUCTIVITY** - This controls 'G' in the conductivity measurement equation mentioned above. 'CW' operation of this adjustment increases 'G'.

**PHASE ERROR** - This lamp comes on when the sensor's electrodes have become damaged or fouled and indicates that the conductivity cell is having difficulty passing current into the solution. Momentary phase errors cause this lamp's output to latch for about one second.

**OVER RANGE** - This lamp comes on when either the conductivity or temperature circuits have voltages within them that are outside of the operating range. Having 'G' settings that are too high usually causes this condition. Momentary over ranges cause this lamp's output to latch for about one second.

The Model 125 back panel contains several connectors. The function of each is fairly obvious except for the circuit outputs that are BNCs. The outer BNC supplies Conductivity  $V_o$ . The inner BNC supplies Temperature  $V_o$ .

## GROUNDING

The Model 125 must always be used with grounded solutions. You have the option of grounding the solution yourself or using the ground supplied by the preamplifier-sensor shaft. The preamplifier-sensor shaft is wired through the cable to JMP 1 on the main bridge's circuit board, near the preamp connector. JMP 1 is accessible by removing the two screws that hold the bridge's cover. The Model 125 is shipped with JMP 1 shorted, connecting the preamp-sensor shaft to circuit ground. If you break this connection, then the preamp-sensor shaft and the sensor float free of circuit ground. This allows you to provide ground to the solution in whatever way that is required. The Model 125 circuit ground must be connected to this ground through the outer side of either BNC on the back panel.

## CALIBRATION

Precision Measurement Engineering (PME) performs temperature calibration when the Model 125 is initially supplied. In this condition, a small piece of tape is placed over the GAIN-TEMPERATURE potentiometer access hole. This is intended to prevent accidental adjustment of this potentiometer creating an uncalibrated condition. The conductivity channel is not calibrated by PME and must be calibrated prior to each use. Note that when a new sensor is installed, the customer must perform both calibrations. Instructions are given below for both temperature and conductivity calibration.

Calibration of the Model 125 is performed in three stages: 'GAIN' setup for proper range, calibration data acquisition, and numerical analysis. Each of these is discussed below:

**GAIN SETUP (G)** - It is important that the instrument does not go out of range during an experiment. To insure that it does not, you must first set both 'GAIN' adjustments at the maximum experimental conditions (Important note: If you wish to continue to use PME's temperature calibration do NOT set the temperature potiotometer. Skip this in the following instructions. See previous paragraph.) The procedure is:

1. Turn both gain pots fully CCW. (Counter ClockWise rotation of at least twenty turns.) The over range lamps for both temperature and conductivity should be dark.
2. Place the sensor in a solution at the lowest temperature expected during the experiment. Adjust

the temperature gain pot CW (ClockWise rotation) until the over range light comes on. Adjust the pot CCW ½ to 1 turn so that the light goes out. (If the light does not come on while adjusting CW more than twenty turns, then leave the pot set fully CW.)

3. Place the sensor in a solution at the highest conductivity expected. Adjust the conductivity gain pot CW until the over range light comes on. Adjust the pot CCW ½ to 1 turn so that the light goes out. (If the light does not come on while adjusting CW more than twenty turns, then leave the post set fully CW.)

**CALIBRATION DATA ACQUISITION-** This process depends on the accuracy that you need in your experiment. Calibration of the Model 125 depends on laboratory standards for temperature and conductivity that you supply. The Model 125 will not measure more accurately than these standards. The procedure for simple calibration using a 0.1 deg C resolution thermometer and one solution of know conductivity is given below.

1. Turn the electronic bridge off. Disconnect the preamplifier-sensor cable. Turn the electronic bridge on. Record the voltage at the temperature output. This voltage is Temperature Voff. It should be approximately -5 volts. Turn the power off and reconnect the cable. Turn the power on.
2. Expose the sensor to two stirred solutions at temperatures approximately bracketing the expected temperature span of the experiment. Record the temperature from the reference thermometer and the voltage at the temperature output. (Additional temperature points within the span may also be recorded.)
3. Hold the sensor with the tip up and gently shake the water from the electrodes. Record the voltage at the conductivity output. This voltage is conductivity Voff. It should be approximately -5 volts.
4. Place the sensor in a solution of known conductivity that is near the highest conductivity expected. (Conductivity can be calculated from density and temperature or salt concentration and temperature or otherwise known.) Record the conductivity and the voltage at the conductivity output. (You should expect to see roughly +4 Volts if conductivity gain is properly adjusted.)

**NUMERICAL ANALYSIS-** A numerical means for connecting voltage outputs from the Model 125 with their corresponding conductivity and temperatures is described below.

1. Use the temperature response equation given in the Specifications section. Include 'G' in the coefficients 'A' and 'B'. Find a simultaneous solution to the equation:

$$\ln (V(T)-V_{\text{off}}) = A + B/T \quad (T \text{ in deg K})$$

for 'A' and 'B' by using the two (optionally more) data points. If your reference measurement of 'T' was absolutely certain, then this equation will give fit errors of about 0.05 deg C over a 20 deg C span.

2. Use the conductivity response equation given in the Specifications section. Solve:

$$V(C) - V_{\text{off}} = G * C$$

for 'G' using the data point measured above. If your measurement of conductivity was absolutely certain, then this equation will give fit errors of less than 1% of the reading over the linear range of conductivity.

## MAINTENANCE

Maintenance of the Model 125 consists of keeping the sensor, preamplifier, and cable assembly free from physical damage, rinsed off, and dry.

From time to time a sensor assembly will need to be replaced. The procedures are (see following diagrams):

1. Remove the outer screw (screw #1) nearest the cable end of the preamplifier. Place it in a cup to prevent loss. It is very small. Refer to Sensor Replacement drawing.
2. Grip the stainless steel shaft and carefully slide the preamplifier cover toward the sensor until all the interior parts are exposed.
3. Clip the six wires (4 - conductivity and 2 - temperature) about 1 cm from each connecting point on the circuit board. Do this so that some of the color from each wire remains connected to the circuit board. This will help when soldering the new wires.

4. Loosen the set-screw (screw #2) in the shaft retaining bracket at the shaft end of the preamplifier circuit board.
5. Withdraw the stainless steel shaft. Slide the preamplifier cover off the stainless steel shaft. Discard the used shaft.
6. Inspect the o-ring at the cable end of the preamplifier for any cuts or debris. If detected, then replace. Next, apply a small amount of silicon oil around the o-ring.
7. Remove the outer screw (screw #3) on the preamplifier cover that holds the nosepiece. Inspect the o-ring for cuts or debris. If detected, then replace. Apply a small amount of silicon oil around the o-ring.
8. The nosepiece has a smaller, inner o-ring. Carefully remove it and inspect for cuts or debris. If detected, then replace. Apply a small amount of silicon oil around the o-ring before replacing it back into the nosepiece.
9. Slide the nosepiece onto the preamplifier cover. Be sure not to cut the outer o-ring with the cover. Install the outer screw (screw #3). Be careful not to over tighten.
10. Insert the o-ring expander into the stainless steel shaft where the wires protrude. Insert the wires through the clear, plastic tubing. The fit is not tight and the expander will fall to the side. See the Sensor Insertion drawing.
11. Apply a small amount of silicon oil to the end of the new stainless steel shaft and to the o-ring expander. **Do not put silicone oil on the clear, plastic tubing since it will be inserted into the shaft.**
12. Insert the wires through the nosepiece and carefully slide the preamplifier cover onto the o-ring expander. While holding the preamplifier cover, position the o-ring expander so that it is flush with the stainless steel shaft. Keep the expander and shaft in a straight, horizontal position. Gently push the preamplifier cover over the expander/shaft section and onto the shaft being careful not to cut the inner o-ring. See the Shaft Insertion drawing.



13. Once the preamplifier cover is on the shaft, remove the o-ring expander and insert the wires through the shaft-retaining bracket. Insert the shaft fully into the bracket and tighten the set-screw (screw #2).
14. One at time, unsolder the old wire on the circuit board and solder the new wire (that has been stripped) onto the connecting point. Be careful not to switch the colors. If you become confused, then refer to the Sensor Connection drawing. If you coat the stripped end of the wire with solder flux, then it becomes much easier to solder.
15. Be sure to inspect the solders and cut any exposed-stripped wire.
16. Carefully slide the preamplifier cover back into place over the circuit board. To insure that the cover does not cut the wires, be sure that they are positioned near the center of the circuit board. Be careful not to cut the o-ring near the cable end as the cover slides over it.
17. Install the outer screw (screw #1). Be careful not to over tighten. Wipe off any remaining silicon oil on the shaft or preamplifier cover.

## REFERENCES

Thermistor Calibration: Thermometrics (1995), **Thermistors Catalog**, Thermometrics, Inc., 808 U.S. Hwy 1, Edison, New Jersey 08817, Phone 1(800) 246-7019. This catalog is free upon request from Thermometrics.

Conductivity for Sodium Chloride Solutions: Chiu, Ying-Check and Fuoss, Raymond M., "Conductance of the Alkali Halides XII. Sodium and Potassium Chloride in Water at 25 C." **Journal of Physical Chemistry**, Vol. 72, # 12, (1968), pp. 4123 - 4129.


Hewitt, G.F., "Tables of the Resistivity of Aqueous Sodium Chloride Solutions", Atomic Energy Research Establishment (Great Britain), **Report R3497**, (1960).

Surdo, Antonio Lo, Alzda, E.M. and Millero, F.J., "The PVT Properties of Concentrated Aqueous Electrolytes. I. Densities and Apparent Molal Volumes of NaCl, Na<sub>2</sub>SO<sub>4</sub>, MgCl<sub>2</sub>, and MgSO<sub>4</sub> solutions from 0.1 to Saturation and from 273.15 to 323.15 K", **Journal of Chemical Thermodynamics**, Vol. 14 (1982), pp. 649 - 662.

Conductivity for Seawater: "Special Issue on the Practical Salinity Scale", **IEEE J. of Oceanic Engineering**, Vol. OE-5, #1 (1980).

MSCTI Instrument: Head, M. J., "The Use of Miniature Four-Electrode Conductivity Probes for High Resolution Measurement of Turbulent Density or

Temperature Variations in Salt-Stratified Water Flow”, **Ph.D. Thesis, University of California, San Diego**, 1983.




INSTALLING THE SHIELD

- 1) LOOSEN THE 4 RETAINING SCREWS
- 2) PLACE SENSOR AGAINST SHIELD AS SHOWN
- 3) ROTATE SENSOR INTO AXIAL ALIGNMENT WITH SHIELD CENTER HOLE
- 4) GENTLY SLIDE SENSOR INTO SHIELD UNTIL WITHIN 1/4" FROM SHIELD END
- 5) GENTLY TIGHTEN ALL 4 RETAINING SCREWS
- 6) WHILE HOLDING SHIELD (NOT SENSOR SHAFT) SLIDE COVER OVER SHIELD
- 7) AVOID SHARP TAPS ON SENSOR - THEY RESONATE THE CONDUCTIVITY SENSOR AND WILL DAMAGE IT.

REV	REVISIONS DESCRIPTION	DATE

PRECISION MEASUREMENT ENGINEERING 1827 HAWK VIEW DRIVE ENCINITAS, CA 92024	
SENSOR SHIELD	4250
DATE	SHEET 1 OF 1

TOLERANCES (unless otherwise noted)
X .XX
.XX .XXX




PRECISION MEASUREMENT ENGINEERING  
 1827 HAWK VIEW DRIVE  
 ENCINITAS, CA 92024

**SENSOR REPLACEMENT**

A | 4215 | SHEET 5 OF 6  
 DATE 20-OCT-91

TOLERANCES (unless otherwise noted)

.x  
 .xx  
 .xxx




PRECISION MEASUREMENT ENGINEERING  
 1827 HAWK VIEW DRIVE  
 RICHMOND, CA 94804

SENSOR INSERTION

DATE 20-OCT-91 SHEET 6 OF 6

TOLERANCES (unless otherwise noted)

.X .XX .XXX +.0000 +.0005



PRECISION MEASUREMENT ENGINEERING  
 1827 HAWK VIEW DRIVE  
 ENCINITAS, CA 92024  
**SENSOR CONNECTION**  
 4215  
 DATE 27-APRIL-89 SHEET 3 OF 3

TOLERANCES (unless otherwise noted)  
 .x  
 .xx  
 .xxx